Extraction of AB from Brain for AB ELISA

This protocol is a method to extract A β from brain for application to A β ELISA

Reagents


- 50mM Tris-buffered saline containing protease inhibitor cocktail, pH 7.6.
- 6M GuHCl solution (6M GuHCl/50mM Tris-HCl, protease inhibitor cocktail, pH 7.6) Store at -30 °C.
- 90% formic acid
- GuHCl diluent (20mM phosphate, 0.4M NaCl, 2mM EDTA, 10% Block Ace, 0.2% BSA, 0.05% NaN3, 0.075% CHAPS, protease inhibitor cocktail, pH 7.0)

Equipment

- Teflon-glass homogenizer (Wheaton 1ml tapered tissue grinder)
- Beckman Optima TL Ultracentrifuge, TLA100.4 rotor
- Motor-driven plastic pestle
- Evaporator

Procedure

Extraction of Aβ from tissues



References

- Saido TC, Iwata N. (2006) Metabolism of amyloid beta peptide and pathogenesis of Alzheimer's disease- Towards presymptomatic diagnosis, prevention and therapy-*Neurosci. Res.* 54(4): 235-253.
- Huang S-M, Mouri A, Kokubo H, Nakajima R, Suemoto T, Higuchi M, Staufenbiel M, Noda Y, Yamaguchi H, Nabeshima T, Saido TC, Iwata N. (2006). Neprilysinsensitive synapse-associated Abeta oligomers impair neuronal plasticity and cognitive function. *J. Biol. Chem.* 281(26): 17941-17951.
- 3 Saito T, Iwata N, Tsubuki S, Takaki Y, Takano J, Huang S-M, Suemoto T, Higuchi M, Saido TC. (2005). Somatostatin regulates brain amyloid beta peptide, Abeta42, via modulation of proteolytic degradation. *Nature Med.* 11(4): 434-439.
- 4 Iwata N, Mizukami H, Shirotani K, Takaki Y, Muramatsu S, Lu B, Gerard NP, Gerard C, Ozawa K, Saido TC. (2004). Presynaptic localization of neprilysin contributes to efficient clearance of amyloid-beta peptide in mouse brain. *J. Neurosci.* 24(4): 991-998.
- 5 Iwata N, Tsubuki S, Takaki Y, Shirotani K, Bao L, Gerard NP, Gerard C, Hama E, Lee H-J, Saido TC. (2001). Metabolic Regulation of Brain Abeta by Neprilysin. Science 292(5521): 1550-1552.